


# Anti-Human CD107b (LAMP-2) FITC

Catalog Number: 11-1078

Also Known As:LAMP2, lysosomal-associated membrane protein 2 RUO: For Research Use Only. Not for use in diagnostic procedures.



Intracellular staining of the Jurkat cell line with Mouse IgG1 K Isotype Control FITC (cat. 11-4714) (open histogram) or Anti-Human CD107b (LAMP-2) FITC (filled histogram). Total cells were used for analysis.

#### **Product Information**

Contents: Anti-Human CD107b (LAMP-2) FITC

REF Catalog Number: 11-1078 Clone: eBioH4B4 (H4B4)

> Concentration: 5 uL (0.25 ug)/test Host/Isotype: Mouse IgG1, kappa

Formulation: aqueous buffer, 0.09% sodium azide, may contain carrier protein/stabilizer

Temperature Limitation: Store at 2-8°C. Do not freeze. Light sensitive material.

LOT Batch Code: Refer to Vial

Use By: Refer to Vial

Contains sodium azide

### Description

The eBioH4B4 monoclonal antibody reacts with human CD107b, also known as lysosomal-associated membrane protein-2 (LAMP-2). CD107b is a highly glycosylated, type I transmembrane protein of approximately 105 kDa. It is expressed intracellularly in lysosomal/endosomal membranes in nearly all cells. It is also expressed on the surface of degranulating T cells (to a lesser extent than CD107a) and activated platelets as well as some cancer cells. In humans, mutations in CD107b results in a lysosomal glycogen storage disorder, known as Danon disease.

#### Applications Reported

This eBioH4B4 (H4B4) antibody has been reported for use in intracellular staining followed by flow cytometric analysis. It has also been reported for use in surface staining in a flow cytometric based degranulation assay.

## **Applications Tested**

This eBioH4B4 (H4B4) antibody has been pre-titrated and tested by intracellular staining and flow cytometric analysis of Jurkat cells. This can be used at 5 µL (0.25 µg) per test. A test is defined as the amount (µg) of antibody that will stain a cell sample in a final volume of 100 µL. Cell number should be determined empirically but can range from 10<sup>5</sup> to 10<sup>8</sup> cells/test.

#### References

Betts MR, Koup RA. Detection of T-cell degranulation: CD107a and b. Methods Cell Biol. 2004;75:497-512. (H4B4, FC, degranulation assay, PubMed)

Grutzkau A, Smorodchenko A, Lippert U, Kirchhof L, Artuc M, Henz BM. LAMP-1 and LAMP-2, but not LAMP-3, are reliable markers for activationinduced secretion of human mast cells. Cytometry A. 2004 Sep;61(1):62-8. (H4B4, FC, PubMed)

Eskelinen EL, Illert AL, Tanaka Y, Schwarzmann G, Blanz J, Von Figura K, Saftig P. Role of LAMP-2 in lysosome biogenesis and autophagy. Mol Biol Cell. 2002 Sep;13(9):3355-68.

Nishino I, Fu J, Tanji K, Yamada T, Shimojo S, Koori T, Mora M, Riggs JE, Oh SJ, Koga Y, Sue CM, Yamamoto A, Murakami N, Shanske S, Byrne E, Bonilla E, Nonaka I, DiMauro S, Hirano M. Primary LAMP-2 deficiency causes X-linked vacuolar cardiomyopathy and myopathy (Danon disease). Nature. 2000;406:906-910. (H4B4, IHC, PubMed)

Carlsson SR, Roth J, Piller F, Fukuda M. Isolation and characterization of human lysosomal membrane glycoproteins, h-lamp-1 and h-lamp-2. Major sialoglycoproteins carrying polylactosaminoglycan. J Biol Chem. 1988 Dec 15;263(35):18911-9.

#### **Related Products**

11-1072 Anti-Mouse CD107b (LAMP-2) FITC (eBioABL-93 (ABL-93)) 11-1079 Anti-Human CD107a (LAMP-1) FITC (eBioH4A3)

Not for further distribution without written consent. Copyright © 2000-2012 eBioscience, Inc.

Tel: 888.999.1371 or 858.642.2058 • Fax: 858.642.2046 • www.eBioscience.com • info@eBioscience.com