

Monoclonal Antibody to Bovine IFN- γ

Antibody:	PAN
Product code:	3115-7
Quantity:	100 tests
Immunogen:	Bovine IFN- γ
Isotype:	Mouse IgG1
Specificity:	Native and recombinant bovine IFN- γ . Cross-reacts with native ovine IFN- γ and native and recombinant equine IFN- γ . In addition, mAb PAN also reacts with porcine IFN- γ .
Contents:	Ready-to-use solution of conjugated antibody in sterile filtered (0.2 μ m) PBS with 0.2% BSA and 0.09% sodium azide. 5 μ l is recommended for staining of 1 million cells in a total volume of 50 μ l.
Purification:	Purified from <i>in vitro</i> cultures by protein G affinity chromatography.
Conjugation:	The PAN antibody has been conjugated to fluorescein isothiocyanate (FITC).
Storage:	Store product at 4-8°C or frozen at -20°C or below. Avoid repeated freezing/thawing.
Applications:	Detection of IFN- γ producing bovine, equine, ovine and porcine cells by flow cytometry.



Note; for research use only.

Mabtech shall not be liable for the use or handling of the product or for consequential, special, indirect or incidental damages therefrom.



MABTECH AB
Box 1233
SE-131 28 Nacka Strand
Sweden
Tel +46 8 716 27 00
Fax +46 8 716 27 01
mabtech@mabtech.com
www.mabtech.com

Mabtech, Inc.
M.E.B. Suite 220
3814 West Street
Mariemont, OH 45227
USA
Toll Free: +1 866 ELI-SPOT
Tel: +1 513 871 4500
Fax: +1 513 871 7353
mabtech.usa@mabtech.com

Mabtech Australia Pty Ltd
44 Gresswell Rd
Macleod, VIC 3085
AUSTRALIA
Tel: +61 3 9459 9630
Fax: +61 3 9455 0084
mabtech.au@mabtech.com

Mabtech AB
Büro Deutschland
GERMANY
Tel: +49 40 4135 7935
Fax: +49 40 4135 7945
mabtech.de@mabtech.com

Mabtech AB
Bureau de liaison France
BP 255, 1300 route des
Crêtes
06905 Sophia Antipolis
FRANCE
Tel: +33 4 92 38 80 70
Fax: +33 4 92 38 80 71
mabtech.fr@mabtech.com