

Product Data Sheet

FITC anti-human CD23

Catalog # / Size: 338505 / 25 tests

338506 / 100 tests

Clone: EBVCS-5 **Isotype:** Mouse $IgG1 \kappa$

Reactivity: Human

Preparation: The antibody was purified by affinity chromatography, and conjugated with

FITC under optimal conditions. The solution is free of unconjugated FITC.

Formulation: Phosphate-buffered solution, pH 7.2, containing 0.09% sodium azide and

0.2% (w/v) BSA (origin USA).

Storage: The antibody solution should be stored undiluted at 4°C and protected from

prolonged exposure to light. Do not freeze.

Applications:

Applications: FC - Quality tested

Recommended Usage: Each lot of this antibody is quality control tested by immunofluorescent

staining with flow cytometric analysis. Test size products are transitioning from 20 µl to 5 µl per test. Please check your vial or your CoA to find the suggested use of this reagent per million cells in 100 µl staining volume or per 100 µl of whole blood. It is recommended that the reagent be titrated for

optimal performance for each application. Read more at www.biolegend.com/testsize regarding the test size change. Application References: 1. Sugden B and Metzenberg S. 1983. J. Virol. 46:800-807.

Description: CD23 is a 45 kD protein, also known as Leu-20, FcɛRII, IgE Fc receptor, BLAST-2, B6, and low affinity IgE receptor. It is a member of the Ig family, expressed on most mature B cells, B cells in follicular mantle (but not in proliferating germinal center cells, follicular dendritic cells, monocytes, eosinophils, Langerhans cells, and a subset of T cells (10-15% of tonsillar T cells). CD23 responds to high levels of IgE by downregulating IgE secretion. In human monocytes, CD23 triggering results in release of pro-inflammatory cytokines including TNF-α, IL-1, IL-6, and GM-CSF. CD23 can be proteolytically cleaved to generate soluble CD23 fragments of various molecular weights. In chronic lymphocytic leukemia, levels of soluble CD23 in the serum can be used as a prognostic marker to identify patients at high risk for disease progression. Alternate splicing of exon 2 can also generate two cell-surface isoforms of CD23 differing by 6 amino acids in their cytoplasmic region. The D3.6 antibody against human CD23 has been shown to be useful for flow cytometry.

Clone

MOPC-21

- Antigen References: 1. Ludin C, et al. 1987. EMBO J. 6:109.
 - Delespesse G, et al. 1992. Immunol. Rev. 125:77.
 Flores-Romo L, et al. 1993. Science 261:1038. 4. Armant M, et al. 1994. J. Exp. Med. 180:1005.

Related Products: Product

FITC Mouse IgG1, κ Isotype Ctrl (FC)

Cell Staining Buffer RBC Lysis Buffer (10X)

Human TruStain FcX™ (Fc Receptor Blocking Solution)

EBVCS-5 FITC Human peripheral blood lymphocytes stained with CD19 APC and

EBVCS-5 FITC

Application

FC, ICC, ICFC FC, ICFC FC, ICC, ICFC



