


An Affymetrix Company

# Anti-Mouse IL-12 p35 Alexa Fluor 647 (To Be Discontinued, Refer to Cat. No. 50-7352)

Catalog Number: 51-7352

Also known as: Interleukin-12 p35

RUO: For Research Use Only. Not for use in diagnostic procedures.



In vitro-cultured mouse dendritic cells unstimulated (left) or stimulated with Mouse IFN gamma Recombinant Protein (cat. 14-8311), LPS, and the Imidazoquinoline, R848 (right), in the presence of Protein Transport Inhibitor Cocktail (cat. 00-4980) stained intracellularly with Anti-Mouse IL-12/IL-23 p40 PE (cat. 12-7123) and 0.06 ug of Anti-Mouse IL-12 p35 Alexa Fluor® 647. Total viable cells were used for analysis.

### **Product Information**


**Contents:** Anti-Mouse IL-12 p35 Alexa Fluor 647 (To Be Discontinued, Refer to Cat. No.

50-7352)

REF Catalog Number: 51-7352

Clone: 4D10p35

Concentration: 0.2 mg/mL Host/Isotype: Rat IgG2a, kappa Formulation: aqueous buffer, 0.09% sodium azide, may contain carrier protein/stabilizer Temperature Limitation: Store at 2-8°C. Do not





### Description

This 4D10p35 monoclonal antibody reacts with mouse p35, a component of two of the four members of the interleukin (IL)-12 family of cytokines. The IL-12 family members are heterodimeric cytokines composed of an alphachain, such as p35, and a beta-chain, such as p40. The alpha-chain proteins are thought to be the limiting factors for secretion of IL-12 family member cytokines, as expression of the alpha-chain is more restricted than that of the beta-chain proteins, and the two subunits may need to be co-expressed in order for cytokine secretion to occur. IL-12 family members are secreted primarily by antigen-presenting cells in response to various TLR agonists and other inflammatory stimuli. The exact types of stimuli present and kinetics of stimulation determine the balance of IL-12 family member production.

The p35 protein was first identified as part of the IL-12 cytokine, which is also known as IL-12p70. Secretion is induced by IFN gamma and agonists of TLR3, 4, or 7 and is inhibited by Th2 cytokines. IL-12 is essential for the differentiation and maintenance of Th1 cells. The other subunit, p40, is expressed in much higher quantities than p35, and may be secreted independently as a monomer or homodimer. Free p40 may act as an antagonist of IL-12 function.

More recently, p35 was also found to dimerize with EBI3 (Epstein-Barr Virus-Induced Gene 3) to form IL-35. The biology of IL-35 is incompletely understood, and it is believed to be unique among the IL-12 family members in its expression and function. It is produced by regulatory T cells, and appears to have immunoregulatory functions. IL-35 has been found to suppress the proliferation and function of effector T cells, and has also been observed to induce a regulatory phenotype on naïve CD4+ cells cultured in its presence.

Testing of 4D10p35 in the context of IL-35 is currently in process. For additional information, please contact technical support.

### **Applications Reported**



An Affymetrix Company

## Anti-Mouse IL-12 p35 Alexa Fluor 647 (To Be Discontinued, Refer to Cat. No. 50-7352)

Catalog Number: 51-7352

Also known as: Interleukin-12 p35

RUO: For Research Use Only. Not for use in diagnostic procedures.

This 4D10p35 antibody has been reported for use in intracellular staining followed by flow cytometric analysis.

### **Applications Tested**

This 4D10p35 antibody has been tested by intracellular staining followed by flow cytometric analysis. This can be used at less than or equal to 0.125  $\mu$ g per test. A test is defined as the amount ( $\mu$ g) of antibody that will stain a cell sample in a final volume of 100  $\mu$ L. Cell number should be determined empirically but can range from 10<sup>5</sup> to 10<sup>8</sup> cells/test. It is recommended that the antibody be carefully titrated for optimal performance in the assay of interest.

#### References

Collison LW, Chaturvedi V, Henderson AL, Giacomin PR, Guy C, Bankoti J, Finkelstein D, Forbes K, Workman CJ, Brown SA, Rehg JE, Jones ML, Ni HT, Artis D, Turk MJ, Vignali DA. IL-35-mediated induction of a potent regulatory T cell population. Nat Immunol. 2010 Dec;11(12):1093-101.

Goriely S., Neurath MF., and Goldman M. How microorganisms tip the balance between interleukin-12 family members. Nat Rev Immunol. 2008 Jan;8(1):81-6

Collison L and Vignali D. Interleukin-35: Odd one out or part of the family? Immunol Rev. 2008 Dec; 226: 248-62.

Yoshimura A, Yoshida H, Miyazaki Y, Kinjyo I, Ishibashi T, Yoshimura T, Takeda A, Hamano S. Two sided roles of IL-27: induction of Th1 differentiation on Naïve CD4+ T cells versus suppression of proinflammatory cytokine production including IL-23-induced IL-17 on activated CD4+ T cells partially through STAT3-dependent mechanism. J Immunol. 2006; 177: 5377-85.

### **Related Products**

00-4980 Protein Transport Inhibitor Cocktail (500X)

12-7123 Anti-Mouse IL-12/IL-23 p40 PE (C17.8)

14-8311 Mouse IFN gamma Recombinant Protein

50-4321 Rat IgG2a K Isotype Control eFluor® 660 (Alexa Fluor® 647 Replacment) (eBR2a)

50-7023 Anti-Mouse IL-23 p19 eFluor® 660 (fc23cpg)

51-7359 Anti-Human IL-12 p35 Alexa Fluor® 647 (To Be Discontinued. Refer to Cat. No. 50-7359) (SNKY35)

88-7120 Mouse IL-12/IL-23 total p40 ELISA Ready-SET-Go!®

88-7121 Mouse IL-12 p70 ELISA Ready-SET-Go!®

88-8823 Intracellular Fixation & Permeabilization Buffer (plus Brefeldin A) (previously named IC Fixation & Permeabilization Buffer)

#### l egal

Alexa Fluor® and Pacific Blue® are registered trademarks of and licensed under patents assigned to Molecular Probes, Inc. for research use only. This product is subject to an agreement between Molecular Probes, Inc. and eBioscience, and the manufacture, use, sale or import of this product may be subject to one or more U.S. patents, pending applications and corresponding foreign equivalents, owned by Molecular Probes, Inc. (a wholly owned subsidiary of Invitrogen Corp). The purchase of this product conveys to the buyer the non-transferable right to use the purchased amount of the product for life science research or as an ASR. The buyer cannot use this product for manufacturing or for any other screening (specifically including use in combination with microarrays or High Content Screening) or testing purpose, other than as an ASR. For information on purchasing a license to this product for purposes other than life science research or use as an ASR, contact Molecular Probes, Inc.