

NuPAGE® MOPS SDS Buffer Kit (for NuPAGE® MES SDS Buffer Kit (for Bis-Tris Gels), contains 1 ea. NP0002, NP0004, NP0005, NP0007) (Novex®)

货号 **NP0060**
包装 **1 kit**

描述

The NuPAGE® MES SDS 缓冲液试剂盒专门为使用 NuPAGE® Novex® Bis-Tris 凝胶分离小到中型分子蛋白而设计。包括以下缓冲液：

- NuPAGE® MES SDS 电泳缓冲液 (20X, 500 ml - NP0002)，用于小到中型分子蛋白分离
- NuPAGE® 样品还原剂 (10X, 250 µl - NP0004)
- NuPAGE® 抗氧化剂 (NP0005)
- NuPAGE® LDS 样品 缓冲液 (4X, 10 ml - NP0007)

预混的 NuPAGE® 缓冲液通过便捷的方式，来确保高品质的，一致性的电泳结果。所有的缓冲液都用高纯度试剂配制，并经过严格的质量控制。缓冲液是浓缩型的，在使用前只需用蒸馏水稀释即可。

仅供研究使用。不适用于动物或人类治疗或诊断用。

详细说明

常用规格

Gel Type:	SDS PAGE Gel
Buffer Type:	Running Buffers
Product Size:	1 kit
Gel Compatibility:	NuPAGE® Bis-Tris Gels
Regulatory Statement:	For Research Use Only. Not for use in diagnostic procedures.

内容和儲存

Contains

- NuPAGE® MES SDS Running Buffer (20X, 500 ml)
- NuPAGE® Sample Reducing Agent (10X, 250 µl) containing 500 mM Dithiothreitol
- NuPAGE® Antioxidant (15 ml)
- NuPAGE® LDS Sample Buffer (10 ml)

Store each at 4°C. Warm NuPAGE® LDS Sample Buffer to 25°C before use.