

Masterflex® Peristaltic Pumps

Product numbers VFP001 & VFP002



Masterflex peristaltic pumps perfectly compliment Vivaflow® cross flow filtration modules, offering the flexibility, durability, and performance expected for today's laboratory applications.

Information

The peristaltic design confines the fluid to the tubing so there is no contact between pump and sample. This means there is no danger of the pump contaminating the sample, and the sample cannot contaminate the pump.

Peristaltic pumps are positive displacement pumps. This means they are self-priming, non-siphoning with no back flow. Unlike gear or piston pumps, peristaltic pumps can be used to pump liquids, gasses or mixed phase fluids.

For ease of cleaning, when changing samples simply change the tubing; with the Easy Load pump head options, this takes a mater of seconds. The painted steel housing can be wiped down with a mild detergent solution.

Product features

- Peristaltic design
- Variable speed
- Sample fluid is confined to the tubing
- Easy to clean
- Positive displacement
- Versatile
- Low maintenance

Specifications

Product number	VFP001	VFP002	
Line voltage (50 – 60 Hz)	180 – 260 V	90 – 130 V	
Drive speed range	20 – 600 rpm	20 – 600 rpm	
Flow range	16 – 480 ml/min (size 16 tubing) 33 – 1,000 ml/min (size 15 tubing)	16 – 480 ml/min (size 16 tubing) 33 – 1,000 ml/min (size 15 tubing)	
Maximum torque	635.5 mN·m	90 oz·in	
Power output	37 W	0.05 hp	
Maximum current	0.9 A	1.5 A	
Dimensions	220 × 180 × 135 mm	811/16 × 71/8 × 55/16 in	
Weight	4.1 kg	9 lb	

Masterflex products at a glance

Ordering information

Product	Pack size	Product number
Masterflex Economy drive variable speed peristaltic pump (240 V)	1	VFP001
Masterflex Economy drive variable speed peristaltic pump (115 V)	1	VFP002
Masterflex Easy Load pump head (size 16)	1	VFA012
Masterflex Easy Load pump head (size 15)	1	VFA013
Size 16 PVC pump tubing and luer fittings (3 m, 3.2×1.6 mm)	1	VFA004
Size 15 PVC pump tubing and luer fittings (3 m, 4.8×2.6 mm)	1	VFA003

Sartorius Stedim Biotech GmbH August-Spindler-Strasse 11 37079 Goettingen, Germany Phone +49.551.308.0 Fax +49.551.308.3289 www.sartorius-stedim.com

USA Toll-Free +1.800.368.7178 UK +44.1372.737159 France +33.442.845600 Italy +39.055.63.40.41 Spain +34.90.2110935 Japan +81.3.3740.5407

Specifications subject to change without notice. Printed and copyrighted by Sartorius Stedim Biotech GmbH. | W Publication No.: SL-2068-e11031 Order No.: 35032-539-57 Ver. 03 | 2011