Technical Data Sheet

PE Mouse Anti-Human Cytokeratin 14, 15, 16 and 19 Set

Product Information

Material Number: 550953 Size: 100 tests

QC Testing: Human Reactivity:

51-38255X Component:

Description: PE Mouse Anti-Human Cytokeratin

Size: 100 tests (1 ea)

Vol. per Test: 20 µl Clone Name: KA4

Isotype: Mouse IgG1, κ

Storage Buffer: Aqueous buffered solution containing BSA and ≤0.09% sodium azide.

Component: 51-35405X-6


Description: PE Mouse IgG1 Isotype Control

100 tests (1 ea) Size: Vol. per Test: 20 µl Clone Name: MOPC-21 Isotype: Mouse IgG1, κ

Storage Buffer: Aqueous buffered solution containing BSA and ≤0.09% sodium azide.

Description

Reacts with cytokeratins 14, 15, 16 and 19 present in the cytoplasm of a wide variety of human simple epithelium. Human epithelium has been shown to contain cytoplasmic filaments of the cytokeratin type. There have been 19 different polypeptides identified and they demonstrate a characteristic pattern for each type of epithelium. This reagent is useful for intracellular staining for flow cytometric analysis, western blot and immunohistochemical staining of both, acetone-fixed frozen tissue sections and formalin-fixed paraffin embedded tissue sections with citrate pretreatment.



Profile of fixed and permeabilized A431 cells with 1% paraformaldehyde and 0.1% Triton X-100 and analyzed by flow cytometry

Preparation and Storage

Store undiluted at 4°C and protected from prolonged exposure to light. Do not freeze.

The monoclonal antibody was purified from tissue culture supernatant or ascites by affinity chromatography.

The antibody was conjugated with R-PE under optimum conditions, and unconjugated antibody and free PE were removed.

Application Notes

Application

Intracellular staining (flow cytometry) Routinely Tested

BD Biosciences

bdbiosciences.com

Europe Japan 32.2.400.98.95 0120.8555.90 Latin America/Caribbean **United States** Asia Pacific 800.268.5430

For country contact information, visit bdbiosciences.com/contact

Conditions: The information disclosed herein is not to be construed as a recommendation to use the above product in violation of any patents. BD Biosciences will not be held responsible for patent infringement or other violations that may occur with the use of our products. Purchase does not include or carry any right to resell or transfer this product either as a stand-alone product or as a component of another product. Any use of this product other than the permitted use without the express written authorization of Becton, Dickinson and Company is stictly prohibited. For Research Use Only. Not for use in diagnostic or therapeutic procedures. Not for resale.

Unless otherwise noted, BD, BD Logo and all other trademarks are property of Becton, Dickinson and Company. © 2014 BD



Recommended Assay Procedure:

Note: This product is routinely tested on fixed and permeabilized A431 cells with 1% paraformaldehyde and 0.1% Triton X-100.

Product Notices

- This antibody has been optimized and preassayed with its matched isotype control to be used at the recommended volume of 20 ul/test. Titration of the reagents or substituting with other (non-matched) isotype control is NOT recommended.
- Caution: Sodium azide yields highly toxic hydrazoic acid under acidic conditions. Dilute azide compounds in running water before discarding to avoid accumulation of potentially explosive deposits in plumbing.
- 3. Source of all serum proteins is from USDA inspected abattoirs located in the United States.
- For fluorochrome spectra and suitable instrument settings, please refer to our Multicolor Flow Cytometry web page at www.bdbiosciences.com/colors.
- Please refer to www.bdbiosciences.com/pharmingen/protocols for technical protocols. 5.

Nagle RB, Lucas DO, McDaniel KM, Clark VA, Schmalzel GM. Paget's cells. New evidence linking mammary and extramammary Paget cells to a common cell phenotype. Am J Clin Pathol. 1985; 83(4):431-438. (Biology)

Nagle RB, Moll R, Weidauer H, Nemetschek H, Franke WW. Different patterns of cytokeratin expression in the normal epithelia of the upper respiratory tract. Differentiation. 1985; 30(2):130-140. (Biology)



bdbiosciences.com

 Canada
 Europe
 Japan

 800.268.5430
 32.2.400.98.95
 0120.8555.90
 Asia Pacific 65.6861.0633 United States Latin America/Caribbean For country contact information, visit bdbiosciences.com/contact

Conditions: The information disclosed herein is not to be construed as a recommendation to use the above product in violation of any patents. BD Biosciences will not be held responsible for patent infringement or other violations that may occur with the use of our products. Purchase does not include or carry any right to resell or transfer this product either as a stand-alone product or as a component of another product. Any use of this product other than the permitted use without the express written authorization of Becton, Dickinson and Company is stictly prohibited.

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Not for resale.

Unless otherwise noted, BD, BD Logo and all other trademarks are property of Becton, Dickinson and Company. © 2014 BD



550953 Rev. 4 Page 2 of 2