Technical Data Sheet

Purified Mouse Anti-6XHIS with Control

Product Information

Material Number: 552565 Size: $50 \, \mu g$


51-9000012 Component: **Description:** 6xHis Antibody Size: 50 μg (1 ea) Clone Name: F24-796 Isotype: Mouse IgG2b, κ

Storage Buffer: Aqueous buffered solution containing ≤0.09% sodium azide.

51-9000129 Component: 6xHis Control Lysate **Description:** Size: 100 μl (1 ea)

Description

The engineering of a short epitope tag to a recombinant protein provides a very convenient method for detecting and isolating the protein of interest. For example, the placement of a series of 6 consecutive histidine residues that are incorporated at either end of the recombinant fusion protein facilitates the purification of that protein. Purification can be easily accomplished by metal chelate affinity chromatography using a metal resin column. The 6xHis monoclonal antibody can be used to detect the epitope-tagged protein of interest.



Western blot analysis. An E coli-produced recombinant 6xHis control lysate was probed with anti-6xHis monoclonal antibody (clone F24-796, Cat. No. 552565) at the following concentrations: 0.5 (lane 1), 0.25 (lane 2), and 0.125 µg/ml (lane 3). 6xHis control is identified as a band of ~31 kDa.

Preparation and Storage

The monoclonal antibody was purified from tissue culture supernatant or ascites by affinity chromatography. Store undiluted at 4°C.

Application Notes

Application

Western blot	Routinely Tested
Immunoprecipitation	Tested During Development

Recommended Assay Procedure:

An E. coli produced recombinant 6xHis control lysate is provided as a positive control. Store lysates at -20°C.

Note: The size of the band observed will vary according to the size of the recombinant protein. In the example shown, the lysate contained a recombinant protein of ~31 kDa and the 6xHis tag does not appreciably change the observed molecular weight.

BD Biosciences

bdbiosciences.com

United States Asia Pacific Latin America/Caribbean Canada Europe 877.232.8995 888.259.0187 32.53.720.550 0120.8555.90 65.6861.0633

For country-specific contact information, visit bdbiosciences.com/how_to_order/

Conditions: The information disclosed herein is not to be construed as a recommendation to use the above product in violation of any patents. BD Biosciences will not be held responsible for patent infringement or other violations that may occur with the use of our products. Purchase does not include or carry any right to resell or transfer this product either as a stand-alone product or as a component of another product. Any use of this product other than the permitted use without the express written authorization of Becton Dickinson and Company is strictly prohibited. For Research Use Only. Not for use in diagnostic or the respect to proceedures. Not for resale.

BD, BD Logo and all other trademarks are the property of Becton, Dickinson and Company. ©2008 BD



Product Notices

- 1. Since applications vary, each investigator should titrate the reagent to obtain optimal results.
- 2. Please refer to www.bdbiosciences.com/pharmingen/protocols for technical protocols.
- Caution: Sodium azide yields highly toxic hydrazoic acid under acidic conditions. Dilute azide compounds in running water before discarding to avoid accumulation of potentially explosive deposits in plumbing.

References

Hochuli E, Dobeli H, Schacher A. New metal chelate adsorbent selective for proteins and peptides containing neighbouring histidine residues. *J Chromatogr.* 1987; 411:177-184.(Biology)

552565 Rev. 3 Page 2 of 2