## **Technical Data Sheet**

# PE Mouse anti-Ezrin (pY353)

#### **Product Information**

558601 **Material Number:** 50 tests Size: 20 µl Vol. per Test: I66-386 Clone:

Phosphorylated Human Ezrin Peptide Immunogen:


Mouse IgG1 κ Isotype: Human Reactivity:

Aqueous buffered solution containing BSA and ≤0.09% sodium azide. Storage Buffer:

#### Description

Ezrin is a member of the ERM (ezrin-radixin-moesin) family of proteins that function as crosslinkers between the actin cytoskeleton and the plasma membrane. Ezrin is a substrate for the tyrosine kinase epidermal (EGF) and hepatocyte growth factor (HGF) receptors. Upon EGF or HGF stimulation, ezrin is phosphorylated on tyrosine 353 (Y353), which is required for its interaction with phosphatidylinositol 3 kinase p85 subunit and activation of Akt.

The I66-386 monoclonal antibody recognizes the phosphorylated Y353 of human ezrin.



Analysis of Ezrin (pY353) in human chronic myelogenous leukemia. Serum-starved K-562 cells (ATCC CCL-243) were either treated with 1 mM Pervanadate (Sigma Cat. No. S6508) at 37°C for 20 minutes (shaded histogram) or untreated (open histogram). The cells were fixed (BD Cytofix™ buffer, Cat. No. 554655) for 10 minutes at 37 °C, then permeabilized (BD™ Phosflow Perm Buffer III, Cat. No. 558050) on ice for at least 30 minutes, and then stained with PE Mouse anti-Ezrin (pY353). Flow cytometry was performed on a BD™ FACSCalibur flow cytometry system.

## **Preparation and Storage**

The monoclonal antibody was purified from tissue culture supernatant or ascites by affinity chromatography.

The antibody was conjugated with R-PE under optimum conditions, and unconjugated antibody and free PE were removed by gel

Store undiluted at 4° C and protected from prolonged exposure to light. Do not freeze.

## **Application Notes**

### Application

Intracellular staining (flow cytometry)	Routinely Tested

## **Suggested Companion Products**

Catalog Number	Name	Size	Clone
554655	Fixation Buffer	100 ml	(none)
558050	Perm Buffer III	125 ml	(none)

## **BD Biosciences**

bdbiosciences.com

United States Europe 877.232.8995 888.259.0187 32.53.720.550 0120.8555.90 65.6861.0633 55.11.5185.9995 For country-specific contact information, visit bdbiosciences.com/how\_to\_order/

Conditions: The information disclosed herein is not to be construed as a recommendation to use the above product in violation Conditions: The information disclosed herein is not to be construed as a recommendation to use the above product in violation drap patents. BD Biosciences will not be held responsible for patent infringement or other violations that may occur with the use of our products. Purchase does not include or carry any right to resell or transfer this product either as a stand-alone product or as a component of another product. Any use of this product other than the permitted use without the express written authorization of Becton Dickinson and Company is strictly prohibited.

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Not for resale.

BD, BD Logo and all other trademarks are the property of Becton, Dickinson and Company. ©2006 BD



### **Product Notices**

- 1. Please refer to www.bdbiosciences.com/pharmingen/protocols for technical protocols.
- 2. This reagent has been pre-diluted for use at the recommended Volume per Test. We typically use 1 X 10e6 cells in a  $100-\mu l$  experimental sample (a test).
- For fluorochrome spectra and suitable instrument settings, please refer to our Fluorochrome Web Page at www.bdbiosciences.com/pharmingen/colors.
- Caution: Sodium azide yields highly toxic hydrazoic acid under acidic conditions. Dilute azide compounds in running water before discarding to avoid accumulation of potentially explosive deposits in plumbing.
- 5. Source of all serum proteins is from USDA inspected abattoirs located in the United States.

#### References

Gautreau A, Poullet P, Louvard D, Arpin M. Ezrin, a plasma membrane-microfilament linker, signals cell survival through the phosphatidylinositol 3-kinase/Akt pathway. *J Biol Chem.* 1999; 96:7300-7305.(Biology)

Krieg J, Hunter T. Identification of the two major epidermal growth factor-induced tyrosine phosphorylation sites in the microvillar core protein ezrin. *J Biol Chem.* 1992; 267(27):19258-19265.(Biology)

558601 Rev. 1 Page 2 of 2