Technical Data Sheet

Biotin Hamster Anti-Rat CD62L

Product Information

559989 **Material Number:**

Alternate Name: L-selectin, LECAM-1

0.1 mg 0.5 mg/mlConcentration: Clone: HRL1

Rat LECAM-1 human IgG1 fusion protein Immunogen:

Armenian Hamster IgG2, λ1 Isotype:

Reactivity: QC Testing: Rat

Storage Buffer: Aqueous buffered solution containing ≤0.09% sodium azide.

Description

The HRL1 antibody reacts with L-selectin (LECAM-1, CD62L), which is detected on a small percentage of thymocytes and on most neutrophils and peripheral lymphocytes. CD62L is a 62-kDa (on neutrophils) or 65-kDa (on lymphocytes) receptor, with lectin-like and epidermal growth factor-like domains, which binds to sialylated oligosaccharide determinants on high endothelial venules (HEV) in peripheral lymph nodes. This member of the selectin adhesion molecule family appears to be required for lymphocyte homing to peripheral lymph nodes and to contribute to neutrophil emigration at inflammatory sites. L-selectin is rapidly shed from lymphocytes and neutrophils upon cell activation. In the mouse, the level of CD62L expression distinguishes naive CD4+T cells from effector/memory T helper cells. HRL1 antibody inhibits the ligand-binding activity of Lselectin in in vitro assays and slows in vivo leukocyte rolling on microvascular endothelium.

Preparation and Storage

The monoclonal antibody was purified from tissue culture supernatant or ascites by affinity chromatography. The antibody was conjugated with biotin under optimum conditions, and unreacted biotin was removed. Store undiluted at 4° C and protected from prolonged exposure to light. Do not freeze.

Application Notes

Application

rippleation			
	Flow cytometry	Routinely Tested	

Recommended Assay Procedure:

For detection of L-selectin on peripheral blood leukocytes, use of EDTA as anti-coagulant is recommended; heparin should not be used.

Suggested Companion Products

Catalog Number	Name	Size	Clone
553963	Biotin Hamster IgG2, λ1 Isotype Control	0.25 mg	Ha4/8
554061	PE Streptavidin	0.5 mg	(none)

Product Notices

- Since applications vary, each investigator should titrate the reagent to obtain optimal results.
- Please refer to www.bdbiosciences.com/pharmingen/protocols for technical protocols.
- Although hamster immunoglobulin isotypes have not been well defined, BD Biosciences Pharmingen has grouped Armenian and Syrian hamster IgG monoclonal antibodies according to their reactivity with a panel of mouse anti-hamster IgG mAbs. A table of the hamster IgG groups, Reactivity of Mouse Anti-Hamster Ig mAbs, may be viewed at http://www.bdbiosciences.com/pharmingen/hamster_chart_11x17.pdf.
- Caution: Sodium azide yields highly toxic hydrazoic acid under acidic conditions. Dilute azide compounds in running water before discarding to avoid accumulation of potentially explosive deposits in plumbing.

BD Biosciences

www.bdbiosciences.com

Asia Pacific Europe 32.53.720.550 0120.8555.90 877.232.8995 888.259.0187 65.6861.0633 55.11.5185.9995 For country-specific contact information, visit www.bdbiosciences.com/how to order/

Conditions: The information disclosed herein is not to be construed as a recommendation to use the above product in violation

Conditions: In einformation disclosed neign is not to be construed as a recommendation to use the above product in violation of any patents. BD Biosciences will not be held responsible for patent infringement or other violations that may occur with the use of our products. Purchase does not include or carry any right to resell or transfer this product either as a stand-alone product or as a component of another product. Any use of this product other than the permitted use without the express written authorization of Becton Dickinson and Company is strictly prohibited.

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Not for resale.

BD, BD Logo and all other trademarks are the property of Becton, Dickinson and Company. @2007 BD



References

Mulligan MS, Miyasaka M, Tamatani T, Jones ML, Ward PA. Requirements for L-selectin in neutrophil-mediated lung injury in rats. J Immunol. 1994;

559989 Rev. 15 Page 2 of 2