

01-6101 10 ml Normal Rabbit Serum

Lot No. See product label

**NORMAL RABBIT SERUM** 



01-6101 10 ml Normal Rabbit Serum

Lot No. See product label

NORMAL RABBIT SERUM

QUANTITY:

10 ml

FORM:

Lyophilized. Reconstitution with 10 ml of distilled water will give a 10 ml whole rabbit serum and 0.1% sodium azide.

STORAGE:

Store at 4°C. After reconstitution, store product in frozen aliquots.

**REMARKS:** 

For Research Use Only. Not for use in diagnostic procedures.

QUANTITY:

10 ml

FORM:

Lyophilized. Reconstitution with 10 ml of distilled water will give a 10 ml whole rabbit serum and 0.1% sodium azide.

STORAGE:

Store at 4°C. After reconstitution, store product in frozen aliquots.

**REMARKS:** 

For Research Use Only. Not for use in diagnostic procedures.

www.invitrogen.com

Invitrogen Corporation • 542 Flynn Rd • Camarillo • CA 93012 • Tel: 800.955.6288 E-mail: techsupport@invitrogen.com

PI016101 (Rev 10/08) DCC-08-1089

Important Licensing Information - These products may be covered by one or more Limited Use Label Licenses (see the Invitrogen Catalog or our website, www.invitrogen.com). By use of these products you accept the terms and conditions of all applicable Limited Use Label Licenses. Unless otherwise indicated, these products are for research use only and are not intended for human or animal diagnostic, therapeutic or commercial use.

www.invitrogen.com

Invitrogen Corporation • 542 Flynn Rd • Camarillo • CA 93012 • Tel: 800.955.6288 E-mail: techsupport@invitrogen.com

PI016101

(Rev 10/08) DCC-08-1089

Important Licensing Information - These products may be covered by one or more Limited Use Label Licenses (see the Invitrogen Catalog or our website, www.invitrogen.com). By use of these products you accept the terms and conditions of all applicable Limited Use Label Licenses. Unless otherwise indicated, these products are for research use only and are not intended for human or animal diagnostic, therapeutic or commercial use.