

Mouse (monoclonal) Anti-Human Interleukin-10 (IL-10)

PRODUCT ANALYSIS SHEET

Catalog Number: AHC9102

Lot Number: See product label

Clone Number: 945A2A5

Quantity/Volume: 0.5 mg/0.5 mL

Form of the Antibody: Purified immunoglobulin in phosphate buffered saline, pH 7.5, with 0.05% sodium azide

(Caution: sodium azide is a poisonous and hazardous substance. Handle with care and

dispose of properly.)

Purification: Purified by Protein A affinity column chromatography.

Immunogen: Purified recombinant human IL-10.

Myeloma/

Fusion Partners: BALB/c mouse spleen cells were fused with SP/20 mouse myeloma cells.

Isotype: IgG1 kappa

Specificity: This antibody recognizes natural and recombinant human IL-10, and cross-reacts with

the Epstein-Barr viral IL-10.

Applications: Suitable for use in immunohistochemistry with cryostat sections. Tonsil tissue can be

used as a positive control. A concentration range of 10 µg-20 µg/mL is recommended for immunohistochemistry. The optimal concentration should be determined for each

specific application.

Dilutions: Further dilutions should be made in medium or buffered solution with carrier protein,

such as PBS with 0.1-1.0% BSA.

Storage: Store at 2-8°C in aliquots for maximal stability.

References: Van Hoffen, E., D. van Wichen, I. Stuij, N. de Jonge, C. de Klopping, J. Lahpor, J. van

den Tweel, F. Gmelig-Meyling, and R. de Weger (1996) In situ expression of cytokines

in human heart allografts. American Journal of Pathology 149:1991-2003.

Hong, Y.S., H. Laks, G.G. Cui, T. Chong, and L.Y. Sen (2002) Localized immunosuppression in the cardiac allograft induced by a new liposome-mediated IL-10

gene therapy. J. Heart Lung Transplantation 21(11):1188-1200.

This product is for research use only. Not for use in diagnostic procedures.

www.invitrogen.com

 $Invitrogen\ Corporation \bullet 542\ Flynn\ Rd \bullet Camarillo \bullet CA\ 93012 \bullet Tel:\ 800.955.6288 \bullet E-mail: \\ \underline{techsupport@invitrogen.com}$

PI AHC9102

(Rev 10/08) DCC-08-1089