

References for Products 17518 and 17519

1. Jarikote DV, Krebs N, Tannert S, Roder B, Seitz O. (2006) Exploring Base-Pair-Specific Optical Properties of the DNA Stain Thiazole Orange. *Chemistry*.
2. Trevithick-Sutton CC, Mikelsons L, Filippenko V, Scaiano JC. (2006) Effect of UVC-Induced Photooxidation of DNA on the Intercalation of Thiazole Orange: A Convenient Reporter for Oxidative Damage. *Photochem Photobiol*.
3. Algar WR, Massey M, Krull UJ. (2006) Fluorescence resonance energy transfer and complex formation between thiazole orange and various dye-DNA conjugates: implications in signaling nucleic acid hybridization. *J Fluoresc*, 16, 555.
4. Karunakaran V, Perez Lustres JL, Zhao L, Ernsting NP, Seitz O. (2006) Large dynamic Stokes shift of DNA intercalation dye Thiazole Orange has contribution from a high-frequency mode. *J Am Chem Soc*, 128, 2954.
5. Skripchenko A, Wagner SJ, Thompson-Montgomery D, Awatefe H. (2006) Thiazole orange, a DNA-binding photosensitizer with flexible structure, can inactivate pathogens in red blood cell suspensions while maintaining red cell storage properties. *Transfusion*, 46, 213.
6. Fechter EJ, Olenyuk B, Dervan PB. (2005) Sequence-specific fluorescence detection of DNA by polyamide-thiazole orange conjugates. *J Am Chem Soc*, 127, 16685.
7. Wang X, Krull UJ. (2005) Synthesis and fluorescence studies of thiazole orange tethered onto oligonucleotide: development of a self-contained DNA biosensor on a fiber optic surface. *Bioorg Med Chem Lett*, 15, 1725.
8. Kohler O, Jarikote DV, Seitz O. (2005) Forced intercalation probes (FIT Probes): thiazole orange as a fluorescent base in peptide nucleic acids for homogeneous single-nucleotide-polymorphism detection. *Chembiochem*, 6, 69.
9. Michimoto T, Okamura T, Suzuki K, Watari T, Kano R, Hasegawa A. (2004) Thiazole orange positive platelets in a dog with Evans' syndrome. *J Vet Med Sci*, 66, 1305.
10. Carreon JR, Mahon KP, Jr., Kelley SO. (2004) Thiazole orange-peptide conjugates: sensitivity of DNA binding to chemical structure. *Org Lett*, 6, 517.
11. Jouin H, Daher W, Khalife J, Ricard I, Puijalon OM, Capron M, Dive D. (2004) Double staining of *Plasmodium falciparum* nucleic acids with hydroethidine and thiazole orange for cell cycle stage analysis by flow cytometry. *Cytometry A*, 57, 34.
12. Kohler O, Seitz O. (2003) Thiazole orange as fluorescent universal base in peptide nucleic acids. *Chem Commun (Camb)*, 2938.
13. Okamura T, Kurashige A, Hanahachi A, Kano R, Tsujimoto H, Hasegawa A. (2003) Thiazole orange-positive platelets in cats with thrombocytopenia induced by cyclophosphamide. *Vet Rec*, 152, 506.
14. Figuerres E, Olszewski M, Kletzel M. (2001) A flow cytometric technique using thiazole orange to detect platelet engraftment following pediatric stem-cell transplants. *Cytotherapy*, 3, 277.
15. Privat E, Melvin T, Merola F, Schweizer G, Prodhomme S, Asseline U, Vigny P. (2002) Fluorescent properties of oligonucleotide-conjugated thiazole orange probes. *Photochem Photobiol*, 75, 201.
16. Hanahachi A, Kano R, Watari T, Hasegawa A. (2002) Thiazole orange-positive platelets in a dog with idiopathic thrombocytopenic purpura. *Vet Rec*, 150, 48.
17. Privat E, Melvin T, Asseline U, Vigny P. (2001) Oligonucleotide-conjugated thiazole orange probes as "light-up" probes for messenger ribonucleic acid molecules in living cells. *Photochem Photobiol*, 74, 532.
18. Privat E, Asseline U. (2001) Synthesis and binding properties of oligo-2'-deoxyribonucleotides covalently linked to a thiazole orange derivative. *Bioconjug Chem*, 12, 757.
19. Boger DL, Tse WC. (2001) Thiazole orange as the fluorescent intercalator in a high resolution fid assay for determining DNA binding affinity and sequence selectivity of small molecules. *Bioorg Med Chem*, 9, 2511.

20. Hanahachi A, Moritomo T, Kano R, Watari T, Tsujimoto H, Hasegawa A. (2001) Thiazole orange-positive platelets in healthy and thrombocytopenic dogs. *Vet Rec*, 149, 122.
21. Toba K, Fuse I, Koike T, Tsuchiyama J, Higuchi W, Niikuni K, Abe T, Yano T, Takahashi M, Shibata A, Aizawa Y. (2000) [Clinical application of reticulated platelet analysis using the matic-modified thiazole orange method and flow cytometry]. *Rinsho Ketsueki*, 41, 1151.
22. Bunkenborg J, Gadjev NI, Deligeorgiev T, Jacobsen JP. (2000) Concerted intercalation and minor groove recognition of DNA by a homodimeric thiazole orange dye. *Bioconjug Chem*, 11, 861.
23. Joutsu-Korhonen L, Sainio S, Riikonen S, Javela K, Teramo K, Kekomaki R. (2000) Detection of reticulated platelets: estimating the degree of fluorescence of platelets stained with thiazole orange. *Eur J Haematol*, 65, 66.
24. Svanvik N, Westman G, Wang D, Kubista M. (2000) Light-up probes: thiazole orange-conjugated peptide nucleic acid for detection of target nucleic acid in homogeneous solution. *Anal Biochem*, 281, 26.
25. Robinson M, MacHin S, Mackie I, Harrison P. (2000) In vivo biotinylation studies: specificity of labelling of reticulated platelets by thiazole orange and mepacrine. *Br J Haematol*, 108, 859.
26. Fujii T, Shimomura T, Fujimoto TT, Kimura A, Fujimura K. (2000) A new approach to detect reticulated platelets stained with thiazole orange in thrombocytopenic patients. *Thromb Res*, 97, 431.
27. Prodhomme S, Demaret JP, Vinogradov S, Asseline U, Morin-Allory L, Vigny P. (1999) A theoretical and experimental study of two thiazole orange derivatives with single- and double-stranded oligonucleotides, polydeoxyribonucleotides and DNA. *J Photochem Photobiol B*, 53, 60.
28. Bunkenborg J, Stidsen MM, Jacobsen JP. (1999) On the sequence selective bis-intercalation of a homodimeric thiazole orange dye in DNA. *Bioconjug Chem*, 10, 824.
29. Bondensgaard K, Jacobsen JP. (1999) Solution structure and energy calculation of bis-intercalation of homodimeric thiazole orange dye derivatives in DNA: effects of modifying the linker. *Bioconjug Chem*, 10, 735.
30. Balduini CL, Noris P, Spedini P, Belletti S, Zambelli A, Da Prada GA. (1999) Relationship between size and thiazole orange fluorescence of platelets in patients undergoing high-dose chemotherapy. *Br J Haematol*, 106, 202.
31. Petersen M, Hamed AA, Pedersen EB, Jacobsen JP. (1999) Bis-intercalation of homodimeric thiazole orange dye derivatives in DNA. *Bioconjug Chem*, 10, 66.
32. Nygren J, Svanvik N, Kubista M. (1998) The interactions between the fluorescent dye thiazole orange and DNA. *Biopolymers*, 46, 39.
33. Staerk D, Hamed AA, Pedersen EB, Jacobsen JP. (1997) Bisintercalation of homodimeric thiazole orange dyes in DNA: effect of modifying the linker. *Bioconjug Chem*, 8, 869.
34. Faridi J, Nielsen KE, Stein PC, Jacobsen JP. (1997) Dynamic bis-intercalation of a homodimeric thiazole orange dye in DNA: evidence of intercalator creeping. *J Biomol Struct Dyn*, 15, 321.
35. Russell KE, Perkins PC, Grindem CB, Walker KM, Sellon DC. (1997) Flow cytometric method for detecting thiazole orange-positive (reticulated) platelets in thrombocytopenic horses. *Am J Vet Res*, 58, 1092.
36. Hansen LF, Jensen LK, Jacobsen JP. (1996) Bis-intercalation of a homodimeric thiazole orange dye in DNA in symmetrical pyrimidine-pyrimidine-purine-purine oligonucleotides. *Nucleic Acids Res*, 24, 859.
37. Dale GL, Friese P, Hynes LA, Burstein SA. (1995) Demonstration that thiazole-orange-positive platelets in the dog are less than 24 hours old. *Blood*, 85, 1822.
38. Jacobsen JP, Pedersen JB, Hansen LF, Wemmer DE. (1995) Site selective bis-intercalation of a homodimeric thiazole orange dye in DNA oligonucleotides. *Nucleic Acids Res*, 23, 753.

39. Wilke WW, Heller MJ, Jakoubova OK, Robinson RA. (1994) Use of thiazole orange homodimer as an alternative to ethidium bromide for DNA detection in agarose gels. *Mod Pathol*, 7, 385.
40. Bonan JL, Rinder HM, Smith BR. (1993) Determination of the percentage of thiazole orange (TO)-positive, "reticulated" platelets using autologous erythrocyte TO fluorescence as an internal standard. *Cytometry*, 14, 690.
41. Serke S, Huhn D. (1993) Improved specificity of determination of immature erythrocytes (reticulocytes) by multiparameter flow-cytometry and thiazole orange using combined staining with monoclonal antibody (anti-glycophorin-A). *Clin Lab Haematol*, 15, 33.
42. Mechetner EB, Sedmak DD, Barth RF. (1991) Heterogeneity of peripheral blood reticulocytes: a flow cytometric analysis with monoclonal antibody HAE9 and thiazole orange. *Am J Hematol*, 38, 61.
43. Abbott DL, McGrath JP. (1991) Evaluation of flow cytometric counting procedure for canine reticulocytes by use of thiazole orange. *Am J Vet Res*, 52, 723.
44. Bowen D, Bentley N, Hoy T, Cavill I. (1991) Comparison of a modified thiazole orange technique with a fully automated analyser for reticulocyte counting. *J Clin Pathol*, 44, 130.
45. Rye HS, Quesada MA, Peck K, Mathies RA, Glazer AN. (1991) High-sensitivity two-color detection of double-stranded DNA with a confocal fluorescence gel scanner using ethidium homodimer and thiazole orange. *Nucleic Acids Res*, 19, 327.
46. Chin-Yee I, Keeney M, Lohmann RC. (1991) Flow cytometric reticulocyte analysis using thiazole orange; clinical experience and technical limitations. *Clin Lab Haematol*, 13, 177.
47. Van Hove L, Goossens W, Van Duppen V, Verwilghen RL. (1990) Reticulocyte count using thiazole orange. A flow cytometry method. *Clin Lab Haematol*, 12, 287.
48. Kienast J, Schmitz G. (1990) Flow cytometric analysis of thiazole orange uptake by platelets: a diagnostic aid in the evaluation of thrombocytopenic disorders. *Blood*, 75, 116.
49. Davis BH, Bigelow NC. (1989) Flow cytometric reticulocyte quantification using thiazole orange provides clinically useful reticulocyte maturity index. *Arch Pathol Lab Med*, 113, 684.
50. Van Bockstaele DR, Peetermans ME. (1989) 1,3'-Diethyl-4,2'-quinolythiacyanine iodide as a "thiazole orange" analogue for nucleic acid staining. *Cytometry*, 10, 214.
51. Carter JM, McSweeney PA, Wakem PJ, Nemet AM. (1989) Counting reticulocytes by flow cytometry: use of thiazole orange. *Clin Lab Haematol*, 11, 267.
52. Makler MT, Lee LG, Recktenwald D. (1987) Thiazole orange: a new dye for Plasmodium species analysis. *Cytometry*, 8, 568.
53. Lee LG, Chen CH, Chiu LA. (1986) Thiazole orange: a new dye for reticulocyte analysis. *Cytometry*, 7, 508.