

References for Product 71

1. Nagatani H, Suzuki S, Fermin DJ, Girault HH, Nakatani K. (2006) Interfacial behavior of sulforhodamine 101 at the polarized water/1,2-dichloroethane interface studied by spectroelectrochemical techniques. *Anal Bioanal Chem*, 386, 633.
2. Nimmerjahn A, Kirchhoff F, Kerr JN, Helmchen F. (2004) Sulforhodamine 101 as a specific marker of astroglia in the neocortex in vivo. *Nat Methods*, 1, 31.
3. Engelhard HH. (1997) Flow cytometric applications of Sulforhodamine 101 as a fluorescent stain for total cellular protein. *Biotech Histochem*, 72, 1.
4. Lefevre C, Kang HC, Haugland RP, Malekzadeh N, Arttamangkul S. (1996) Texas Res-X and rhodamine Red-X, new derivatives of sulforhodamine 101 and lissamine rhodamine B with improved labeling and fluorescence properties. *Bioconjug Chem*, 7, 482.
5. Ulrich W. (1992) Simultaneous measurement of DAPI-sulforhodamine 101 stained nuclear DNA and protein in higher plants by flow cytometry. *Biotech Histochem*, 67, 73.