

References for Product 66

1. Wang ZQ, Liao J, Diwu Z. (2005) N-DEVD-N'-morpholinecarbonyl-rhodamine 110: novel caspase-3 fluorogenic substrates for cell-based apoptosis assay. *Bioorg Med Chem Lett*, 15, 2335.
2. Boonacker E, Elferink S, Bardai A, Fleischer B, Van Noorden CJ. (2003) Fluorogenic substrate [Ala-Pro]2-cresyl violet but not Ala-Pro-rhodamine 110 is cleaved specifically by DPPiV activity: a study in living Jurkat cells and CD26/DPPiV-transfected Jurkat cells. *J Histochem Cytochem*, 51, 959.
3. Zhang X, Wang H, Fu NN, Zhang HS. (2003) A fluorescence quenching method for the determination of nitrite with Rhodamine 110. *Spectrochim Acta A Mol Biomol Spectrosc*, 59, 1667.
4. Cai SX, Zhang HZ, Guastella J, Drewe J, Yang W, Weber E. (2001) Design and synthesis of rhodamine 110 derivative and caspase-3 substrate for enzyme and cell-based fluorescent assay. *Bioorg Med Chem Lett*, 11, 39.
5. Hug H, Los M, Hirt W, Debatin KM. (1999) Rhodamine 110-linked amino acids and peptides as substrates to measure caspase activity upon apoptosis induction in intact cells. *Biochemistry*, 38, 13906.
6. Jeannot V, Salmon JM, Deumie M, Viallet P. (1997) Intracellular accumulation of rhodamine 110 in single living cells. *J Histochem Cytochem*, 45, 403.
7. Ganesh S, Klingel S, Kahle H, Valet G. (1995) Flow cytometric determination of aminopeptidase activities in viable cells using fluorogenic rhodamine 110 substrates. *Cytometry*, 20, 334.
8. Klingel S, Rothe G, Kellermann W, Valet G. (1994) Flow cytometric determination of cysteine and serine proteinase activities in living cells with rhodamine 110 substrates. *Methods Cell Biol*, 41, 449.
9. Morliere P, Mangel WF, Santus R, Huppe G, Reyftmann JP, Kohen E. (1987) Interaction of tetrapyrrolic rings with rhodamine 110 and 123 and with rhodamine 110 derivatives bearing a peptidic side chain. *Biochem Biophys Res Commun*, 146, 107.